



mesohyal™ MELILOT

OEDEMATOUS CELLULITE WITH A VASCULAR COMPONENT.

mesohyal MELILOT is an anticellulite solution combining melilot extract and troxerutin with a hyaluronic acid base.

Its active substances act synergistically for the purpose of improving microflow and reducing fluid retention:

- **Melilot extract** enhances lymph draining and activates return venous flow. By reducing capillary permeability, it reduces oedema and the feeling of swelling.
- **Troxerutin** minimises capillary reactivity, performing an anti-redness action. It also increases collagen synthesis in the vascular walls and has anti-free radical action.

Based on the diagnosis of the patient's aesthetic disorder, **mesohyal MELILOT** is the ideal complement to prepare cocktails with other products from the **mesohyal®** range:

- **mesohyal HYALURONIC**: with skin moisturising and rejuvenating effect.
- **mesohyal NCTC 109**: with intensive biorevitalising effect.
- **mesohyal VITAMINC**: with antioxidant effect.
- **mesohyal DMAE**: with firming effect.
- **mesohyal ORGANIC SILICON**: to regenerate and restructure the skin tissue.
- **mesohyal BIOTIN**: to reactivate cell metabolism.
- **mesohyal CARNITINE**: to increase energy output.
- **mesohyal ARTICHOKE**: to activate lipolytic processes.

For cocktails with **mesohyal HYALURONIC**, load this first because of its high viscosity. A homogeneous mixture will thus be obtained.

The **TREATMENT** is applied over 20 to 30 minutes, and the patient can return immediately to their activities of daily living. The therapeutic protocol will be based on the patient's degree of cellulite and skin condition, requiring 8 to 12 sessions at intervals of at least 1 week.

Although the reducing action of **mesohyal MELILOT** is immediate and long-lasting, it does not stop the effects of bad habits regarding figure, so regular maintenance treatments will be required every 6/8 months to preserve the skin's healthy appearance. Based on the result obtained, this period may be shortened or prolonged.



mesohyal MELILOT is supplied in 2 ml ampoules. Each ampoule is for single use and for a single patient.



Specifications:

Composition:	Melilot extract (0.3%), troxerutin (3%), and mineral salts Sodium hyaluronate: · Biotechnological origin / <i>non cross-linked</i> · Concentration: 2,5 mg/ml
Indications of use:	Oedematous cellulite Fluid retention
Treatment areas:	Sweet clover and troxerutin extract is especially indicated to treat the areas most susceptible to suffering circulatory dysfunction: anterior side of the thighs, buttocks, abdomen, back, and arms.
Presentation:	Class I 2 ml ampoule (monosilicate glass) Package containing 20 ampoules
Storage conditions:	Keep under 25°C Do not expose to high or low temperatures Avoid direct exposure to light
Injection techniques:	Sterile syringes and needles for mesotherapy techniques / manual or multipuncture gun
Injection depth:	Superficial dermis